## Version 1.0 (2016-11-15)

```
TBseq release
```

The TBseq/ folder contains the following files:

- 1. LICENCE.pdf (v1.0)
- 2. README.pdf (v1.0)
- 3. VERSION.pdf (v1.0)
- 4. TBseq.pl (v1.0)

The TBseq/lib/ folder contains the following packages:

- 1. TBreads (v1.0)
- 2. TBbwa (v1.0)
- 3. TBmerge (v1.0)
- 4. TBrefine (v1.0)
- 5. TBpile (v1.0)
- 6. TBtable (v1.0)
- 7. TBvariants (v1.0)
- 8. TBjoin (v1.0)
- (12)
- 9. TBamend (v1.0)
- 10. TBstrains (v1.0)
- 11. TBgroups (v1.0)

The TBseq/opt/ folder contains the following programs:

- 1. sam-tools (v1.3.1)
- 2. sambamba (v0.6.5)
- 3. IGV-tools (v2.3.88)
- 4. picard-tools (v2.7.1)
- 5. BWA (v0.7.15)
- 6. GATK (v3.6)

The TBseq/var/ folder contains the following sub directories and files:

- 1. cat/
  - a. MTB Gene Categories.txt (v1.0)
- 2. ref/
  - a. M.\_abscessus\_CIP-104536T\_2014-02-03.fasta
    - ${\tt M.\_abscessus\_CIP-104536T\_2014-02-03\_genes.txt}$
  - b. M. chimaera DSM44623 2016-01-28.fasta
    - M.\_chimaera\_DSM44623\_2016-01-28\_genes.txt
  - c. M.\_fortuitum\_CT6\_2016-01-08.fasta
    - M.\_fortuitum\_CT6\_2016-01-08\_genes.txt
  - d. M.\_tuberculosis\_H37Rv\_2015-11-13.fasta
    - M.\_tuberculosis\_H37Rv\_2015-11-13\_genes.txt
- 3. res/
  - a. Base\_Calibration\_List.vcf (v1.0)
  - b. Base\_Calibration\_List.vcf.idx (v1.0)
  - c. MTB Extended Resistance Mediating.txt (v1.0)
  - d. MTB\_Resistance\_Mediating.txt (v1.0)

The tbstest/ folder contains the following files:

- 1. 0002-01\_lib0002\_miseq\_r0001\_151bp\_R2.fastq.gz
- 2. 0002-01\_lib0002\_miseq\_r0001\_151bp\_R1.fastq.gz
- 3. 0001-01-lib0001\_S01\_L004\_R2\_001.fastq.gz
- 4. 0001-01-lib0001 S01 L004 R1 001.fastq.gz

- 5. 0001-01-lib0001\_S01\_L003\_R2\_001.fastq.gz
- 6. 0001-01-lib0001\_S01\_L003\_R1\_001.fastq.gz
- 7. 0001-01-lib0001\_S01\_L002\_R2\_001.fastq.gz
- 8. 0001-01-lib0001\_S01\_L002\_R1\_001.fastq.gz
- 9. 0001-01-lib0001\_S01\_L001\_R2\_001.fastq.gz
- 10.0001-01-lib0001\_S01\_L001\_R1\_001.fastq.gz